GOVERNMENT OF MEGHALAYA HEALTH & FAMILY WELFARE DEPARTMENT

Orders by the Governor NOTIFICATION

Dated Shillong, the 30th July, 2020.

No.Health.141/2020/88:-The Protocol for Rapid Antigen Testing for COVID–19 in Meghalaya is issued for information and compliance of all concerned.

Protocol for Rapid Antigen Testing for COVID-19 in Meghalaya

Detailed instruction

- ❖ Each district to set up adequate number of Antigen Testing Centers (ATC), preferably near containment zones.
- ❖ Temporary arrangement can be made with following condition.
 - ✓ The ATC room should be air-conditioned to ensure that the preferred temperature of the kits is maintained between 2° to 30° C.
 - ✓ ATC should have water for hand washing, electricity & internet.
 - ✓ ATC should have 2-4 rooms or dedicated space divided by partition for arranging sample collection, testing & data entry.
 - ✓ Each room should preferably have one counter and a maximum of 2 counters for the same.
 - ✓ Each ATC should have one qualified Medical doctor for supervision.
 - ✓ Each counter should have Lab Technicians /Rashtriya Bal Swasthya Karyakram (RBSK) Medical Officers, one Data entry operator (DEO) / staff nurse (SN) for entry, dedicated staff to avoid over crowding.

Essentials required for Rapid antigen test:-

- PPE is mandatory for all personnel involved in Sample testing
- Disinfectants/ Hand sanitizer
- Soap and water (for hand washing)
- * Refrigerator/Vaccine carrier box for storing antigen kits
- 1% sodium hypochlorite for cleaning work surfaces
- ❖ Appropriate biomedical waste management bags for disposal of infectious wastes
- VTM and swabs (for collection of samples for RT-PCR for symptomatic negative suspected case

Other instructions

- The rapid antigen test should be performed as per the manufacturer's instructions.
 The video link can be accessed at
 - https://www.youtube.com/watch?v=syXd7kgLSN8
 - https://www.youtube.com/watch?v=mBdaOHJWxI4&feature=youtu.be


- Sample collection should be done by trained personnel wearing proper PPE.
- Only nasopharyngeal swab should be collected from each person. No other sample should be collected.
- ❖ For discarding of the kits and PPE, proper biomedical waste management bags should be available at the ATC.


Recording & reporting in ICMR portal

- > Details of the patient tested must be filled up in the Specimen referral form (SRF).
- > Current version of the SRF available at icmr.gov.in must be used.
- > Each district will be given five login credentials to enter the antigen test details.
- ➤ All the results should be uploaded on the ICMR portal: "https://cvstatus.icmr.gov.in/login.php" on a real time basis.
- Results should be certified by MO
- > No results should go unrecorded
- ➤ In antigen testing, once MO confirms the result, the test shall not be repeated under any circumstances unless advised by MO.
- > SRF for RT-PCR should be duly signed by MO

Interpretation of results, and further action

- ❖ The first reading should be taken at 15 minutes after inoculation of the sample in the cassette.
- For all the samples, the control line should be visible after 15 minutes. If the control line is not visible, then the test should be considered as invalid.
- For positive samples (Red/Pink line-faint), the test line should be visible along with the control line.
- ❖ In case of negative results at 15 minutes, keep the cassette till 30 minutes for final reading.
- ❖ For symptomatic people who are negative by rapid antigen testing, nasal swab/throat swab should be collected in VTM for real-time PCR testing at the nearby lab. The VTM should be sent under cold conditions to an ICMR approved testing laboratory.





Sd/(Sampath Kumar, lAS)
Commissioner & Secretary to the Government of Meghalaya
Health & Family Welfare Department.